

UBA52 Rabbit mAb

货号: B29251

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC ICC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	15kDa
实测分子量	10kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,PC-3,U-937,Mouse liver,Mouse testis,Mouse kidney,Rat liver,Rat testis
细胞定位	cytoplasm,cytosol,cytosolic ribosome,endoplasmic reticulum,extracellular exosome,extracellular space,m itochondrial outer membrane,nucleoplasm,nucleus,plasma membrane
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human UBA52.
序列	MQIFVKLTGTGKTITLEVEPSDTIENVKAKIQDKEGIPPDQQRLLIFAGKQLEDGRTLSDYNIQKESTLHLVLRRLRGGIIEPSLRQLAQKYNCDKMICRKCY

靶点信息

研究背景	Ubiquitin is a highly conserved nuclear and cytoplasmic protein that has a major role in targeting cellular proteins for degradation by the 26S proteasome. It is also involved in the maintenance of chromatin structure, the regulation of gene expression, and the stress response. Ubiquitin is synthesized as a precursor protein consisting of either polyubiquitin chains or a single ubiquitin moiety fused to an unrelated protein. This gene encodes a fusion protein consisting of ubiquitin at the N terminus and ribosomal protein L40 at the C terminus, a C-terminal extension protein (CEP). Multiple processed pseudogenes derived from this gene are present in the genome.
基因ID	7311
基因名	UBA52
Swiss	P62987
别名	UBA52;CEP52;HUBCEP52;L40;RPL40

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn