

IBA1 Rabbit mAb

货号: B29261

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	10kDa/14kDa/16kDa
实测分子量	17kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	THP-1
细胞定位	Cell projection,Cytoplasm,Cytoplasmic side,Peripheral membrane protein,cytoskeleton,ruffle membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human IBA1.
序列	MSQTRDLQG GKAFGLLKAQQEERLDEINKQFLDDPKYSSDEDLP SKLEGFKEKYMEFDLNGNGDIDIMSLKRMLEKLGVP KTHLELKKLIGEVS SSGSETFSYPDFLRMMLGKRSAILKMILMYEEKAREKEKPTGPPAKKAISELP

靶点信息

研究背景	This gene encodes a protein that binds actin and calcium. This gene is induced by cytokines and interferon and may promote macrophage activation and growth of vascular smooth muscle cells and T-lymphocytes. Polymorphisms in this gene may be associated with systemic sclerosis. Alternative splicing results in multiple transcript variants, but the full-length and coding nature of some of these variants is not certain.
------	--

基因ID	199
基因名	AIF1
Swiss	P55008
别名	AIF1;AIF-1;IBA1;IRT-1;IRT1;lba1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn