

Bmi1 Rabbit mAb

货号: B29266

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	36kDa
实测分子量	40kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,A-549,K-562,Mouse heart,Mouse thymus,Rat thymus,Rat testis
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human Bmi1.
序列	LKDYYTLMDIAYIYTWRRNGPLPLKYRVRPTCKRMKISHQRDGLTNAGELESDSGSDKANS PAGGIPSTSSCLPSPSTPVQSPHPQFPHISSTMNGTSNSP

靶点信息

研究背景	This gene encodes a ring finger protein that is major component of the polycomb group complex 1 (PRC1). This complex functions through chromatin remodeling as an essential epigenetic repressor of multiple r egulatory genes involved in embryonic development and self-renewal in somatic stem cells. This protein also plays a central role in DNA damage repair. This gene is an oncogene and aberrant expression is asso ciated with numerous cancers and is associated with resistance to certain chemotherapies. A pseudogen e of this gene is found on chromosome X. Read-through transcription also exists between this gene and t he upstream COMM domain containing 3 (COMMD3) gene.
基因ID	648
基因名	
Swiss	P35226
别名	BMI1;FLVI2/BMI1;PCGF4;RNF51;flvi-2/bmi-1;Bmi1

产品验证



Western blot analysis of Bmi1 expressed in Rat thymus,HeLa,A-549 using Bmi1 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn