

UNG Rabbit mAb

货号: B29277

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	33kDa/34kDa
实测分子量	32kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse skeletal muscle,Mouse heart,Rat skeletal muscle
细胞定位	Mitochondrion,Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human UNG.
序列	MIGQKTLYSFFSPSPARKRHAPSPEPAVQGTGVAGVP EESGDAAAIPAKKAPAGQEEP GTPPSSPLSAEQLDRIQRNKAAA LLRLAARNVPVGFGESWKK

靶点信息

研究背景	This gene encodes one of several uracil-DNA glycosylases. One important function of uracil-DNA glycosylases is to prevent mutagenesis by eliminating uracil from DNA molecules by cleaving the N-glycosylic bond and initiating the base-excision repair (BER) pathway. Uracil bases occur from cytosine deamination or misincorporation of dUMP residues. Alternative promoter usage and splicing of this gene leads to two different isoforms: the mitochondrial UNG1 and the nuclear UNG2. The UNG2 term was used as a previous symbol for the CCNO gene (GeneID 10309), which has been confused with this gene, in the literature and some databases.
基因ID	7374
基因名	UNG
Swiss	P13051
别名	UNG;DGU;HIGM4;HIGM5;UDG;UNG1;UNG15;UNG2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn