

Phospho-Lyn (Y396) Rabbit mAb

货号: **B29284**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	56kDa/58kDa
实测分子量	56kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,HeLa,K-562,Mouse lung,Rat lung
细胞定位	Cell membrane,Cytoplasm,Golgi apparatus,Lipid-anchor,Membrane,Nucleus,perinuclear region
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Phospho-Lyn (Y396).
序列	MGCIKSKGKDSLDDGVDLKTQPVRNTERTIYVRDPTSNKQQRVPVPESQLLPGQRFQTKDPPEEQGDIVVALYPYDGIHPD DLSFKKGKEMKVLEEHGEWWKAKSLLTKKEGFIPSNYVAKLNTLETEEWFFKDITRKDAERQLLAPGNSAGAFLIRESLTK GSFSLSVRDFDPVHGDVIKHYKIRSLDNNGGYISPRITFPCISDMIKHYQKQADGLCRRLEKACISPKPQKPWDKDAW

靶点信息

研究背景	This gene encodes a tyrosine protein kinase, which maybe involved in the regulation of mast cell degranulation, and erythroid differentiation. Alternatively spliced transcript variants encoding different isoforms h ave been found for this gene.
基因ID	4067

基因名	LYN
Swiss	P07948
别名	LYN;JTK8;p53Lyn;p56Lyn

产品验证



Western blot analysis of Phospho-Lyn (Y396) expressed in Mouse lung,U-87MG,HeLa using Phospho-Lyn (Y396) Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn