

TMEM173 Rabbit mAb

货号: B29303

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	
实测分子量	37kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	THP-1
细胞定位	autophagosome,cytosol,endosome,mitochondrial outer membrane,nucleoplasm,perinuclear region of cytoplasm,peroxisome,plasma membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human TMEM173.
序列	LISGM

靶点信息

研究背景	This gene encodes a five transmembrane protein that functions as a major regulator of the innate immune response to viral and bacterial infections. The encoded protein is a pattern recognition receptor that detects cytosolic nucleic acids and transmits signals that activate type I interferon responses. The encoded protein has also been shown to play a role in apoptotic signaling by associating with type II major histocompatibility complex. Mutations in this gene are the cause of infantile-onset STING-associated vasculopathy. Alternate splicing results in multiple transcript variants.
------	--

基因ID	340061
基因名	TMEM173
Swiss	Q86WV6
别名	TMEM173;ERIS;MITA;MPYS;NET23;SAVI;STING;hMITA;hSTING

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn