

PHF8 Rabbit mAb

货号: B29311

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IF/ICC IP
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	140kDa
实测分子量	140kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,C2C12,Mouse thymus
细胞定位	nuclear membrane,nucleolus,nucleoplasm,nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human PHF8.
序列	TVSNSPASQRTPGKRPIKRPAYWRTESEEEEEENASLDEQDSLGA CFKDAEYIYPSLESDDDDPALKSRPKKKKNSDDAPWSPKARVTPTLPKQDRPVREGT

靶点信息

研究背景	The protein encoded by this gene is a histone lysine demethylase that preferentially acts on histones in the monomethyl or dimethyl states. The encoded protein requires Fe(2+) ion, 2-oxoglutarate, and oxygen for its catalytic activity. Defects in this gene are a cause of mental retardation syndromic X-linked Siderius type (MRXSSD). Four transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, May 2010]
------	--

基因ID	23133
基因名	PHF8
Swiss	Q9UPP1
别名	JHDM1F; KDM7B; MRXSSD; ZNF422

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn