

GOLPH2 Rabbit mAb

货号: B29312

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	44kDa/45kDa
实测分子量	45kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,Rat heart
细胞定位	Golgi apparatus,Single-pass type II membrane protein,cis-Golgi network membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human GOLPH2.
序列	TGERLIRVLQDQLKTLQRNYGRLQQDVLQFQKNQTNLERKFSYDLSQCINQMKEVKEQCEERIEEVTKKGNEAVASRDLS ENNDQRQQLQALSEPQPRLQAAGLPHTVEVPQGKGNVLGNSKSQTPAPSSSEVVLDSKRQVEKEETNEIQVVNEEPQRDRLL

靶点信息

研究背景	The Golgi complex plays a key role in the sorting and modification of proteins exported from the endoplasmic reticulum. The protein encoded by this gene is a type II Golgi transmembrane protein. It processes proteins synthesized in the rough endoplasmic reticulum and assists in the transport of protein cargo through the Golgi apparatus. The expression of this gene has been observed to be upregulated in response to viral infection. Alternatively spliced transcript variants encoding the same protein have been described for this gene.
------	---

基因ID	51280
基因名	GOLM1
Swiss	Q8NBJ4
别名	GOLM1;C9orf155;GOLPH2;GP73;HEL46;PSEC0257;bA379P1.3

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn