

NR1D1 Rabbit mAb

货号: B29348

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	66kDa
实测分子量	58kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,Mouse liver,Rat brain,Rat lung
细胞定位	cytoplasm,dendrite,dendritic spine,nuclear body,nucleoplasm,nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human NR1D1.
序列	SPEPTVEDVISQVARAHREIFTYAHDKLGSSPGNFNANHASGSPATTPHRWENQGCPPAPNDNNTLAAQRHNEALNGLRQ

靶点信息

研究背景	This gene encodes a transcription factor that is a member of the nuclear receptor subfamily 1. The encoded protein is a ligand-sensitive transcription factor that negatively regulates the expression of core clock proteins. In particular this protein represses the circadian clock transcription factor aryl hydrocarbon receptor nuclear translocator-like protein 1 (ARNTL). This protein may also be involved in regulating genes that function in metabolic, inflammatory and cardiovascular processes.
------	--

基因ID	9572
基因名	NR1D1
Swiss	P20393
别名	EAR1; hRev; THRA1; THRAL; ear-1; REVERBA; REVERBalpha

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn