

Fascin Rabbit mAb

货号: B29378

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC IP FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50 FC: 1:20 - 1:50
理论分子量	54kDa
实测分子量	54kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,HeLa,NCI-H460,HepG2,SH-SY5Y,Mouse brain,Mouse spinal cord
细胞定位	Cell junction,Cell projection,Cytoplasm,cytoskeleton,cytosol,filopodium,invadopodium,microvillus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human Fascin.
序列	FIGCRKVTGTL DANRSSYDVFQLEFNDGAYNIKDSTGKYWTVGSDSAVTSSGDTPVDFEFEC DYNKVAIKVGGRYLKGD HAGVLKASAETVDPASLWEY

靶点信息

研究背景	This gene encodes a member of the fascin family of actin-binding proteins. Fascin proteins organize F-actin into parallel bundles, and are required for the formation of actin-based cellular protrusions. The encoded protein plays a critical role in cell migration, motility, adhesion and cellular interactions. Expression of this gene is known to be regulated by several microRNAs, and overexpression of this gene may play a role in the metastasis of multiple types of cancer by increasing cell motility. Expression of this gene is also a marker for Reed-Sternberg cells in Hodgkin's lymphoma. A pseudogene of this gene is located on the long arm of chromosome 15.
基因ID	6624
基因名	FSCN1
Swiss	Q16658
别名	FSCN1;FAN1;HSN;SNL;p55;fascin

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn