

Apg12(Atg12) Rabbit mAb

货号: B29404

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC ICC IP FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 ICC: 1:50 - 1:200 IP: 1:20 - 1:50 FC: 1:20 - 1:50
理论分子量	22,55kDa
实测分子量	22,55kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	PC-3,Mouse brain,Mouse stomach
细胞定位	Cytoplasm,Peripheral membrane protein,Preautophagosomal structure membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Apg12(Atg12).
序列	AVGDTPIMKTKKWAVERTRTIQGLIDFIKKFLKLVASEQLFIYVNQSFAPSPDQEVGTLYECFGSDGKLVLHYCKSQAWG

靶点信息

研究背景	Autophagy is a process of bulk protein degradation in which cytoplasmic components, including organelle s, are enclosed in double-membrane structures called autophagosomes and delivered to lysosomes or va cuoles for degradation. ATG12 is the human homolog of a yeast protein involved in autophagy (Mizushim a et al., 1998 [PubMed 9852036]).[supplied by OMIM, Mar 2008]
------	--

基因ID	9140
基因名	ATG12
Swiss	O94817
别名	ATG12;APG12;APG12L;FBR93;HAPG12;Apg12

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn