

Dlx5 Rabbit mAb

货号: B29416

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	20kDa/31kDa
实测分子量	32kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,293T,U-251MG
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Dlx5.
序列	EMPPEHSPSSSDPMACNSPQSPAVWEPQGSSRSLSHHPHAHPPTSNQSPASSYLENSASWY TSAASSINSHLPPPGSLQ HPLALASGTLY

靶点信息

研究背景	This gene encodes a member of a homeobox transcription factor gene family similiar to the Drosophila di stal-less gene. The encoded protein may play a role in bone development and fracture healing. Mutation i n this gene, which is located in a tail-to-tail configuration with another member of the family on the long a rm of chromosome 7, may be associated with split-hand/split-foot malformation.
------	---

基因ID	1749
基因名	DLX5
Swiss	P56178
别名	DLX5;SHFM1D

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn