

KAT13A / SRC1 Rabbit mAb

货号: B29511

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	152kDa/156kDa
实测分子量	157kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	NIH/3T3,K-562
细胞定位	Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human KAT13A / SRC1.
序列	SVTGQVQGNSSIKLELDASKKKESKDHQLRLRYLLDKDEKDLRSTPNLSLDDVKVKVEKKEQMDPCNTNPTPMTKPTPEEI KLEAQSQFTADLDQFDQLPTLEKAAQLPGLCETDRMDGAVTSVTIKSEILPASLQSATARPTSLNRLPELELEAIDNQFG QPGTGDQIPWTNNTVTAINQSKSEDQCISSQLDELLCPPTTVEGRNDEKALLEQLVSFLSGKDETELAELDRALGIDKLVQ GGGLDVLSERFPPQATP

靶点信息

研究背景	The protein encoded by this gene acts as a transcriptional coactivator for steroid and nuclear hormone receptors. It is a member of the p160/steroid receptor coactivator (SRC) family and like other family members has histone acetyltransferase activity and contains a nuclear localization signal, as well as bHLH and PAS domains. The product of this gene binds nuclear receptors directly and stimulates the transcriptional activities in a hormone-dependent fashion. Alternatively spliced transcript variants encoding different isoforms have been identified.
基因ID	8648
基因名	NCOA1
Swiss	Q15788
别名	NCOA1;F-SRC-1;KAT13A;RIP160;SRC1;bHLHe42;bHLHe74

产品验证



Western blot analysis of KAT13A / SRC1 expressed in NIH/3T3,K-562 using KAT13A / SRC1 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn