

gamma Catenin Rabbit mAb

货号: B29533

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	81kDa
实测分子量	82kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,293T,Mouse lung
细胞定位	Cell junction,Cytoplasm,Membrane,Peripheral membrane protein,adherens junction,cytoskeleton,desmosome
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human gamma Catenin.
序列	MEVMNLMEQPIKVTEWQQTYTYDSGIHSGANTCVPSVSSKGIMEEDEACGRQYTLKKTTTYTQGVPPSQGDLEYQMSTT ARAKRVREAMCPGVSGEDSSLLLATQVEGQATNLQRLAEPSQLLSAIVHLINYQDDAELATRALPELTKLLNDEDPVVVTK AAMIVNQLSKKEASRRALMGSPQLVAAVVRTMQNTSDLDTARCTTSILHNLSHHREGLLAIFKSGGIPALVRMLSSPVESVL FYAITTLHNLLLYQEGAKMAVRLADGLQKMVPLLNNPKFLAITTDCLQLLAYGNQ

靶点信息

研究背景	This gene encodes a major cytoplasmic protein which is the only known constituent common to submembranous plaques of both desmosomes and intermediate junctions. This protein forms distinct complexes with cadherins and desmosomal cadherins and is a member of the catenin family since it contains a distinct repeating amino acid motif called the armadillo repeat. Mutation in this gene has been associated with Naxos disease. Alternative splicing occurs in this gene; however, not all transcripts have been fully described.
基因ID	3728
基因名	JUP
Swiss	P14923
别名	JUP;ARVD12;CTNNG;DP3;DP111;PDGB;PKGB

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn