

PAI1 Rabbit mAb

货号: B29553

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IF/ICC IP FC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50 FC: 1:20 - 1:50
理论分子量	43kDa/45kDa
实测分子量	45kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat placenta
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human PAI1.
序列	SEVERARFIINDWVKTHTKGMISNLLGKGAVDQLTRLVLVNALYFNGQWKTPFPDSSTHRRLFHKSDGSTVSVPMMAQT NKFNYTEFTTPDGHHYDILELP

靶点信息

研究背景	This gene encodes a member of the serine proteinase inhibitor (serpin) superfamily. This member is the p rincipal inhibitor of tissue plasminogen activator (tPA) and urokinase (uPA), and hence is an inhibitor of fib rinolysis. Defects in this gene are the cause of plasminogen activator inhibitor-1 deficiency (PAI-1 deficien cy), and high concentrations of the gene product are associated with thrombophilia. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
------	--

基因ID	5054
基因名	SERPINE1
Swiss	P05121
别名	PAI;PAI-1;PAI1;PLANH1;SERPINE1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn