

# c-Jun Rabbit mAb

货号: B29569

产品信息

|       |                                                                                                             |
|-------|-------------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                             |
| 宿主    | Rabbit                                                                                                      |
| 克隆性   | Monoclonal                                                                                                  |
| 预测反应  |                                                                                                             |
| 应用    | WB IHC ICC IP                                                                                               |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IHC:</b> 1:50 - 1:200<br><b>ICC:</b> 1:50 - 1:200<br><b>IP:</b> 1:20 - 1:50 |
| 理论分子量 | 35kDa                                                                                                       |
| 实测分子量 | 43kDa                                                                                                       |
| 形式    | Liquid                                                                                                      |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.               |
| 偶联物   | Unconjugated                                                                                                |
| 阳性对照  | HeLa,HT-1080,NIH/3T3                                                                                        |
| 细胞定位  | Nucleus                                                                                                     |
| 纯化    | Affinity purification                                                                                       |

抗原信息


|      |                                                                                                                                                 |
|------|-------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein corresponding to Human c-Jun.                                                                                        |
| 序列   | TTPTPTQFLCPKNVTDEQEGFAEGFVRALAEHLSQNTLPSVTSAAQPVNGAGMVAPAVASVAGGSGSGGFSASLHSEPP<br>VYANLSNFPNGALSSGGGAPSYGAAGLAFPAQPQQQQPPHLPQQMPVQHPRQLALKEEPQ |

靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                                                       |
|------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene is the putative transforming gene of avian sarcoma virus 17. It encodes a protein which is highl y similar to the viral protein, and which interacts directly with specific target DNA sequences to regulate gene expression. This gene is intronless and is mapped to 1p32-p31, a chromosomal region involved in b oth translocations and deletions in human malignancies. |
|------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |                      |
|-------|----------------------|
| 基因ID  | 3725                 |
| 基因名   | JUN                  |
| Swiss | P05412               |
| 别名    | AP-1;AP1;c-Jun;JUN;c |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)