

GART Rabbit mAb

货号: B29571

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	46kDa/107kDa
实测分子量	107kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,HepG2,HeLa,K-562,Mouse liver
细胞定位	cytosol,extracellular exosome
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human GART.
序列	MAARVLIIGSGGREHTLAWKLAQSHHVKQVLVAPGNAGTACSEKISNTAISISDHTALAQFCKEKKIEFVVVGPEAPLAAGIV GNLRSAGVQCFGPTAEAAQLESSKRFAKEFMDRHGIPTAQWKAFTKPEEACSFILSADFPALVVKASGLAAGKGVIVAKSK EEACKAVQEIMQEKAFGAAGETIVIEELLDGEEVSCLCFTDGKTVAPMPPAQDHKRLLEGDGGPNTGGMGAYCPAPQVSN DLLLLKIKDRTLQRTVDGMQQEGTPYTGILYAGIMLTKNPKVLEFNCRFGDPECQVILPLLKSDLYEVIQSTLDGL

靶点信息

研究背景	The protein encoded by this gene is a trifunctional polypeptide. It has phosphoribosylglycinamide formyltr ansferase, phosphoribosylglycinamide synthetase, phosphoribosylaminoimidazole synthetase activity whi ch is required for de novo purine biosynthesis. This enzyme is highly conserved in vertebrates. Alternativ e splicing of this gene results in two transcript variants encoding different isoforms.
------	---

基因ID	2618
基因名	GART
Swiss	P22102
别名	GART;AIRS;GARS;GARTF;PAIS;PGFT;PRGS

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn