

GCDFP 15 Rabbit mAb

货号: B29575

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IF/ICC IP
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	16kDa
实测分子量	17kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human GCDFP 15.
序列	QDNTRKIIIKNFDIPKSVRPNDEVTAVLAVQTELKECMVVKTYLISSIPLQGAFNYKYTACLDDNPKTFYWDFYTNRTVQIA AVVDVIRELGICPDDAAVIPIKNNRFYTIEILKVE

靶点信息

研究背景	Enables IgG binding activity; aspartic-type endopeptidase activity; and identical protein binding activity. I nvolved in several processes, including detection of chemical stimulus involved in sensory perception of b itter taste; negative regulation of T cell apoptotic process; and proteolysis. Located in extracellular space and nucleus.
基因ID	5304

基因名	PIP
Swiss	P12273
别名	PIP;GCDFP-15;GCDFP15;GPIP4

产品验证



Western blot analysis of GCDFP 15 expressed in HUVEC using GCDFP 15 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn