

SNF2H Rabbit mAb

货号: B29599

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	121kDa
实测分子量	122kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	PC-3,H460,293T,HepG2,SH-SY5Y,Mouse spleen,Mouse brain,Mouse testis
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human SNF2H.
序列	ELFAHFIQPAAQKTPTSPLKMKPGRPRIKKDEKQNLLSVG DYRHR RTEQEDEELLTESSKATNVCTRFEDSPSYVKWGKL RDYQVRGLNWLISLYENGINGILADEMGLGKTLQTISLLGYMKHYRNIPGPHMVLVPKSTLHNWMSEFKRWVPTLRVCL GDKEQRAAFVRDVLLPGEWDVCVTSYEMLIKESVFKFNWRYLVIDEAHRIKNEKSKLSEIVREFKTTNRLLLTGTP LQNN LHELWSLLNFFLPDVFN SADD FDSWFD TNNCLGDQKLVERLHMVLRPFLLRRIKADV

靶点信息

研究背景	The protein encoded by this gene is a member of the SWI/SNF family of proteins. Members of this family have helicase and ATPase activities and are thought to regulate transcription of certain genes by altering the chromatin structure around those genes. The protein encoded by this gene is a component of the chromatin remodeling and spacing factor RSF, a facilitator of the transcription of class II genes by RNA polymerase II. The encoded protein is similar in sequence to the Drosophila ISWI chromatin remodeling protein.
------	--

基因ID	8467
基因名	SMARCA5
Swiss	O60264
别名	SMARCA5;ISWI;SNF2H;WCRF135;hISWI;hSNF2H

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn