

RhoGAP Rabbit mAb

货号: B29604

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	26kDa/27kDa/172kDa
实测分子量	190kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Jurkat,293T,HepG2,PC-12,Mouse brain,Rat brain
细胞定位	Cytoplasm,Membrane,Peripheral membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human RhoGAP.
序列	EGLYRVSGNKTDQDNIQKQFDQDHNINLVSM EVT VNAVAGALKAFFADLPDPLIPYSLHPELLEAAKIPDKTERLHALKEIV KKFHPVNYDVFRYVITHLNRVSQQHKINLMTADNLSICFWPTLMRPDFENREFLSTTKIHQSVVETFIQQCQFFFYNGEIVET TNIVAPPPPSNPGQLVEPMVPLQLPPLQPQLIQPQLQTDPLGII

靶点信息

研究背景	Rho GTPase activating protein 5 negatively regulates RHO GTPases, a family which may mediate cytoskel eton changes by stimulating the hydrolysis of bound GTP. Two transcript variants encoding different isofo rms have been found for this gene.
基因ID	394

基因名	ARHGAP5
Swiss	Q13017
别名	ARHGAP5;GFI2;RhoGAP5;p190-B;p190BRhoGAP

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn