

GATA4 Rabbit mAb

货号: B29608

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	44kDa
实测分子量	54kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SKOV3,Mouse testis,Rat stomach
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human GATA4.
序列	SSNATTSSSEEMRPIKTEPGLSSHGHSSSVSQTFSVSAMSGHGPSIHPVLSALKLSPQGYASPVSQSPQTSSKQDSWNS LVLADSHGDIITA

靶点信息

研究背景	This gene encodes a member of the GATA family of zinc-finger transcription factors. Members of this family recognize the GATA motif which is present in the promoters of many genes. This protein is thought to regulate genes involved in embryogenesis and in myocardial differentiation and function, and is necessary for normal testicular development. Mutations in this gene have been associated with cardiac septal defects. Additionally, alterations in gene expression have been associated with several cancer types. Alternative splicing results in multiple transcript variants.
------	--

基因ID	2626
基因名	GATA4
Swiss	P43694
别名	GATA4;ASD2;TACHD;TOF;VSD1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn