

GC / VDBP Rabbit mAb

货号: B29609

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	39kDa/52kDa/55kDa
实测分子量	53kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Human serum,Mouse serum,THP-1
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human GC / VDBP.
序列	AMDVVFCTYFMPAAQLPELPDVELPTNKDVCDPGNTKVM DKYTFELSRRTHLPEVFLSKVLEPTL KSLGEC CDVEDSTTCF NAKG PLLKKELSSFIDKGQELCADYSENTFTEYKKKLAERLKAKLPDATPTELAKLVNKHSD FASNCCSINSPLYCDSEIDA ELKNIL

靶点信息

研究背景	The protein encoded by this gene belongs to the albumin gene family. It is a multifunctional protein found in plasma, ascitic fluid, cerebrospinal fluid and on the surface of many cell types. It binds to vitamin D and its plasma metabolites and transports them to target tissues. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
------	--

基因ID	2638
基因名	GC
Swiss	P02774
别名	GC;DBP;DBP/GC;GRD3;Gc-MAF;GcMAF;HEL-S-51;VDBG;VDBP

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn