

TAGLN Rabbit mAb

货号: B29611

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	22kDa
实测分子量	23kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse uterus,Rat uterus
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human TAGLN.
序列	MANKGPSYGMSREVQSKIEKKYDEELEERLVEWIIVQCGPDVGRPDGRGLGFQVWLKNGVILSKLVNSLYPDGSKPVKVP ENPPSMVFKQMEQVAQFLKAAEDYGVIKTDMFQTVDLFEGKDMAAVQRTLALGSLAVTKNDGHYRGDPNWFMKKAQ EHKREFTESQLQEGKHVIGLQMGSNRGASQAGMTGYGRPRQIIS

靶点信息

研究背景	The protein encoded by this gene is a transformation and shape-change sensitive actin cross-linking/gelling protein found in fibroblasts and smooth muscle. Its expression is down-regulated in many cell lines, and this down-regulation may be an early and sensitive marker for the onset of transformation. A functional role of this protein is unclear. Two transcript variants encoding the same protein have been found for this gene.
------	--

基因ID	6876
基因名	TAGLN
Swiss	Q01995
别名	TAGLN;SM22;SMCC;TAGLN1;WS3-10

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn