

PGP9.5 Rabbit mAb

货号: B29627

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	25kDa
实测分子量	24kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,DU145,SH-SY5Y,A-549,Mouse brain,Rat brain
细胞定位	Cytoplasm,Endoplasmic reticulum membrane,Lipid-anchor
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human PGP9.5.
序列	MLKLPMEINPEMLNKVLSRLGVAGQWRFVDVLGLEEESLGSPAPACALLLLFPLTAQHENFRKKQIEELKGQEVSPKVYF MKQTIGNSCGTIGLIHAVANNQDKLGFEDGSVLKQFLSETEKMSPEDRAKCFEKNEAIQAAHDAVAQEGQCRVDDKVN HFILFNNVDGHLIELDGRMPFPVNHGASSEDTLKDAAKVCREFTEREQGEVRFSAVALC

靶点信息

研究背景	The protein encoded by this gene belongs to the peptidase C12 family. This enzyme is a thiol protease th at hydrolyzes a peptide bond at the C-terminal glycine of ubiquitin. This gene is specifically expressed in t he neurons and in cells of the diffuse neuroendocrine system. Mutations in this gene may be associated with Parkinson disease.[provided by RefSeq, Sep 2009]
------	---

基因ID	7345
基因名	UCHL1
Swiss	P09936
别名	HEL-117;HEL-S-53;NDGOA;PARK5;PGP 9.5;PGP9.5;PGP95;SPG79;Uch-L1;UCHL1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn