

AKAP 95 Rabbit mAb

货号: B29659

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	76kDa
实测分子量	95kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BT-474,HL-60,SKOV3,HeLa
细胞定位	Nucleus matrix
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human AKAP 95.
序列	SGEPAEDEGPTDTAEAGSDPQAEQLLEEQVPCGTAHEKGVPKARSEAAEAGNGAETMAAEAESAQTRVAPAPAAADAEV EQTDAESKDAVPTE

靶点信息

研究背景	This gene encodes a member of the A-kinase anchor protein family. A-kinase anchor proteins are scaffold proteins that contain a binding domain for the RI/RII subunit of protein kinase A (PKA) and recruit PKA and other signaling molecules to specific subcellular locations. This gene encodes a nuclear A-kinase anchor p rotein that binds to the RII alpha subunit of PKA and may play a role in chromosome condensation during mitosis by targeting PKA and the condensin complex to chromatin. A pseudogene of this gene is located on the short arm of chromosome 9.
------	--

基因ID	10270
基因名	AKAP8
Swiss	O43823
别名	AKAP8;AKAP 95;AKAP-8;AKAP-95;AKAP95

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn