

Bag1 Rabbit mAb

货号: B29667

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	25kDa/31kDa/34kDa/38kDa
实测分子量	33,46,50kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Jurkat,Mouse lung,Mouse testis,Rat testis
细胞定位	Cytoplasm,Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Bag1.
序列	WSEETQSEETQGEEMNRSQEVTRDEESTRSEEVTR EEMAAAGLTVTVTHSNEKHDLHVTSQQGSSEPVVQDLAQVV EEVIGVPQSFQKLIFKGKSLKEM

靶点信息

研究背景	The oncogene BCL2 is a membrane protein that blocks a step in a pathway leading to apoptosis or programmed cell death. The protein encoded by this gene binds to BCL2 and is referred to as BCL2-associated thanogene. It enhances the anti-apoptotic effects of BCL2 and represents a link between growth factor receptors and anti-apoptotic mechanisms. Multiple protein isoforms are encoded by this mRNA through the use of a non-AUG (CUG) initiation codon, and three alternative downstream AUG initiation codons. A related pseudogene has been defined on chromosome X.
------	---

基因ID	573
基因名	BAG1
Swiss	Q99933
别名	BAG1;BAG-1;HAP;RAP46

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn