

# CEND1 Rabbit mAb

货号: B29670

产品信息

|       |                                                                                                                |
|-------|----------------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                                |
| 宿主    | Rabbit                                                                                                         |
| 克隆性   | Monoclonal                                                                                                     |
| 预测反应  |                                                                                                                |
| 应用    | WB IHC IF/ICC IP                                                                                               |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IHC:</b> 1:50 - 1:200<br><b>IF/ICC:</b> 1:50 - 1:200<br><b>IP:</b> 1:20 - 1:50 |
| 理论分子量 | 15kDa                                                                                                          |
| 实测分子量 | 22kDa                                                                                                          |
| 形式    | Liquid                                                                                                         |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.                  |
| 偶联物   | Unconjugated                                                                                                   |
| 阳性对照  | U-251MG,U-87MG,Mouse brain,Rat brain                                                                           |
| 细胞定位  |                                                                                                                |
| 纯化    | Affinity purification                                                                                          |

抗原信息


|      |                                                                                                                                 |
|------|---------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein corresponding to Human CEND1.                                                                        |
| 序列   | MESRGKSASSPKPDTKVPQVTTEAKVPPAADGKAPLTKPSKKEAPAEKQQPPAAPTTAPAKKTSAKADPALLNNHSNLKP<br>APTVPSSPDATPEPKGPGDGAEEDEAASGGPGGRGPWSCENFN |

靶点信息

|      |                                                                                                                                                                                                                                                                        |
|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | The protein encoded by this gene is a neuron-specific protein. The similar protein in pig enhances neurobl<br>astoma cell differentiation in vitro and may be involved in neuronal differentiation in vivo. Multiple pseud<br>ogenes have been reported for this gene. |
| 基因ID | 51286                                                                                                                                                                                                                                                                  |

|       |            |
|-------|------------|
| 基因名   | CEND1      |
| Swiss | Q8N111     |
| 别名    | BM88;CEND1 |

产品验证



Western blot analysis of CEND1 expressed in Mouse brain,Rat brain,U-251MG using CEND1 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)