

ERAP1 Rabbit mAb

货号: B29696

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	100kDa
实测分子量	100kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,Jurkat,Mouse lung,Mouse liver,Rat liver
细胞定位	cytoplasm,cytosol,endoplasmic reticulum,endoplasmic reticulum lumen,extracellular exosome,extracellular region,extracellular space,plasma membrane
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human ERAP1.
序列	NPVGYPLAWQFLRKNWNKLVQKFELGSSSIAHVMVGTNNQFSTRTRLEEVKGFFSSLKENGSQLRCVQQTETIEENIGWMDKNFDKIRVWLQSEKLERM

靶点信息

研究背景	The protein encoded by this gene is an aminopeptidase involved in trimming HLA class I-binding precursors so that they can be presented on MHC class I molecules. The encoded protein acts as a monomer or as a heterodimer with ERAP2. This protein may also be involved in blood pressure regulation by inactivation of angiotensin II. Three transcript variants encoding two different isoforms have been found for this gene. [provided by RefSeq, Oct 2010]
------	---

基因ID	51752
基因名	ERAP1
Swiss	Q9NZ08
别名	A-LAP; ALAP; APPILS; ARTS-1; ARTS1; ERAAP; ERAAP1; PILS-AP; PILSAP

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn