

Desmin Rabbit mAb

货号: B29741

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	53kDa
实测分子量	53kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	C2C12,RD,Mouse lung,Mouse heart,Rat heart
细胞定位	Cell membrane,Cytoplasm,sarcolemma
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Desmin.
序列	EEEIRHLKDEMARHLREYQDLLNVKMALDVEIATYRKLLLEGESRINLPIQTYALSALNFRETSPEQRGSEVHTKKTVMIKTIETRDGEVVSEATQQQHEVL

靶点信息


研究背景	This gene encodes a muscle-specific class III intermediate filament. Homopolymers of this protein form a stable intracytoplasmic filamentous network connecting myofibrils to each other and to the plasma membrane. Mutations in this gene are associated with desmin-related myopathy, a familial cardiac and skeletal myopathy (CSM), and with distal myopathies. [provided by RefSeq, Jul 2008]
------	---

基因ID	1674
基因名	DES
Swiss	P17661
别名	CSM1; CSM2; LGMD2R

产品验证



Immunohistochemical analysis of paraffin-embedded human prostate, using Desmin Antibody.



Western blot analysis of Desmin expressed in Mouse lung,Rat heart,C2C12 using Desmin Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn