

Drebrin Rabbit mAb

货号: B29758

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IF/ICC IP
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	71kDa/76kDa
实测分子量	120kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BT-474,HeLa,SKOV3
细胞定位	Cell junction,Cell projection,Cytoplasm,cell cortex,growth cone
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Drebrin.
序列	EQIERALDEVTSQQPPPLPPPPPAQETQEPSPILDSEETRAAAPQAWAGPMEEPPQAQAPPRGPGSPAEDLMFMESAEQA VLAAPVEPATADATEIHDAADTIETDTATADTTVANNVPPAATSLIDLWPGNGEGASTLQGEPRAPTPPSGTEVTLAEVPLL DEVAPEPLLPAGEGCATLLNFDELPEPPATFCDPEEVEGEPLAAPQTPTLPSALEELEQEQEPEPHLLTNGETTQKEGTQASE GYFSQSQEEFAQSEELCAKAPPPVFYNKPEIDITCWDADPVPEEEEGFEGGD

靶点信息

研究背景	The protein encoded by this gene is a cytoplasmic actin-binding protein thought to play a role in the process of neuronal growth. It is a member of the drebrin family of proteins that are developmentally regulated in the brain. A decrease in the amount of this protein in the brain has been implicated as a possible contributing factor in the pathogenesis of memory disturbance in Alzheimer's disease. At least two alternative splice variants encoding different protein isoforms have been described for this gene.
------	---

基因ID	1627
基因名	DBN1
Swiss	Q16643
别名	DBN1;D0S117E;drebrin

产品验证



Western blot analysis of Drebrin expressed in BT-474,HeLa,SKOV3 using Drebrin Rabbit mAb at 1:1000 . Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blo cking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn