

# RanBP9 Rabbit mAb

货号: B29786

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IF/ICC
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IF/ICC:</b> 1:50 - 1:200
理论分子量	91kDa
实测分子量	91kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	C2C12,RD,C6,Mouse testis,Rat testis
细胞定位	Cytoplasm,Nucleus,cytosol
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human RanBP9.
序列	SGPPAPPGLAAGPGPAGGAPTPALVAGSSAAAPFPHGDSALNEQEKEQLQRRCLKRLYPVDEQETPLPRSWSPKDKFSYIGLSQNNLRVHYKGHGKTPKDAA

靶点信息

研究背景	This gene encodes a protein that binds RAN, a small GTP binding protein belonging to the RAS superfamily that is essential for the translocation of RNA and proteins through the nuclear pore complex. The protein encoded by this gene has also been shown to interact with several other proteins, including met proto-oncogene, homeodomain interacting protein kinase 2, androgen receptor, and cyclin-dependent kinase 11.
基因ID	10048

基因名	RANBP9
Swiss	Q96S59
别名	RANBP9; BPM-L; BPM90; RANBPM; RanBP7; ran-binding protein 9

产品验证



Western blot analysis of RanBP9 expressed in Mouse testis,Rat testis,C2C12 using RanBP9 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)