

# ENTPD5 Rabbit mAb

货号: B29793

## 产品信息

|       |                                                                                               |
|-------|-----------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                               |
| 宿主    | Rabbit                                                                                        |
| 克隆性   | Monoclonal                                                                                    |
| 预测反应  |                                                                                               |
| 应用    | WB IHC IP                                                                                     |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IHC:</b> 1:50 - 1:200<br><b>IP:</b> 1:20 - 1:50               |
| 理论分子量 | 47kDa                                                                                         |
| 实测分子量 | 47kDa                                                                                         |
| 形式    | Liquid                                                                                        |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                  |
| 阳性对照  | Mouse liver,Rat liver,Rat kidney                                                              |
| 细胞定位  | endoplasmic reticulum,extracellular region,extracellular space                                |
| 纯化    | Affinity purification                                                                         |

## 抗原信息


|      |                                                                                                           |
|------|-----------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein corresponding to Human ENTPD5.                                                 |
| 序列   | FEPCYAEVLRVVRGKLHQPEEVQRGSFYAFSYYYDRAVDTDMIDYEKGGILKVEDFERKAREVCDNLENFTSGSPFLCMDL<br>SYITALLKDGFGFADSTVLQ |

## 靶点信息

|      |                                                                                                                                                                                                                                                                                              |
|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | The protein encoded by this gene is similar to E-type nucleotidases (NTPases)/ecto-ATPase/apyrases. NTP ases, such as CD39, mediate catabolism of extracellular nucleotides. ENTPD5 contains 4 apyrase-conserv ed regions which is characteristic of NTPases. [provided by RefSeq, Jan 2009] |
| 基因ID | 957                                                                                                                                                                                                                                                                                          |

|       |                         |
|-------|-------------------------|
| 基因名   | ENTPD5                  |
| Swiss | O75356                  |
| 别名    | CD39L4; NTPDase-5; PCPH |

产品验证



实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)