

NCF1 Rabbit mAb

货号: B29808

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP
推荐浓度	WB: 1:500 - 1:2000 IP: 1:20 - 1:50
理论分子量	30kDa/44kDa
实测分子量	45kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Raji,Daudi,Jurkat(Negative control)
细胞定位	Cytoplasm,Cytoplasmic side,Membrane,Peripheral membrane protein,cytosol
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human NCF1.
序列	MGDTFIRHIALLGFEKRFVPSQHYVYMFLVKWQDLSEKVVYRRFTEIYEFHKTLEMFPIEAGAINPENRIIPHLPAKWFDG QRAAENRQGTLT EYCSLMSLPTKISRCPHLLDFFKVRPDDLKLPDNDQTKKPETYLMPKDGKSTATDITGP IILQTYRAIAN YEKTS GSEMALSTG DVVEVVEKSES GWWFCQMKAKRGWIPASFLEPLDSPDETEDPEPNYAGEPYVAIKAYTAVEGDEV SLLEGEAVEVIHKLLD GWWVIRKDDVTGYFPSMYLQKSGQDVSQAQRQIKRGAPPRRSSIRNAHSIHQSRKRKLSQDAYR RNSVRFLQRRRQARPGPQSPGSPLEERQTQRSKPQPAVPPRPSADLILNRCSESTKRKLASAV

靶点信息

研究背景	The protein encoded by this gene is a 47 kDa cytosolic subunit of neutrophil NADPH oxidase. This oxidase is a multicomponent enzyme that is activated to produce superoxide anion. Mutations in this gene have been associated with chronic granulomatous disease.
------	--

基因ID	653361
基因名	NCF1
Swiss	P14598
别名	NCF1;NCF1A;NOXO2;SH3PXD1A;p47phox

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn