

HMGA1 Rabbit mAb

货号: B29810

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	17kDa
实测分子量	17kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,HepG2,BxPC-3,Raji,Mouse testis
细胞定位	Chromosome,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human HMGA1.
序列	MSESSSKSSQPLASKQEKDGTEKRGRGRPRKQPPVSPGTALVGSQKEPSEVPTPKRPRGRPKGSKNKGAATRKTITTP GRKPRGRPKKLEKEEEEGISQESSEEEQ

靶点信息

研究背景	This gene encodes a chromatin-associated protein involved in the regulation of gene transcription, integr ation of retroviruses into chromosomes, and the metastatic progression of cancer cells. The encoded prot ein preferentially binds to the minor groove of AT-rich regions in double-stranded DNA. Multiple transcript variants encoding different isoforms have been found for this gene. Pseudogenes of this gene have been identified on multiple chromosomes. [provided by RefSeq, Jan 2016]
------	--

基因ID	3159
基因名	HMGA1
Swiss	P17096
别名	HMG-R; HMGA1A; HMGIY

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn