

# DCR2 Rabbit mAb

货号: B29819

## 产品信息

|       |                                                                                               |
|-------|-----------------------------------------------------------------------------------------------|
| 反应    | Human                                                                                         |
| 宿主    | Rabbit                                                                                        |
| 克隆性   | Monoclonal                                                                                    |
| 预测反应  |                                                                                               |
| 应用    | WB IHC IF/ICC                                                                                 |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IHC:</b> 1:50 - 1:200<br><b>IF/ICC:</b> 1:50 - 1:200          |
| 理论分子量 | 41kDa                                                                                         |
| 实测分子量 | 36kDa                                                                                         |
| 形式    | Liquid                                                                                        |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                  |
| 阳性对照  | THP-1                                                                                         |
| 细胞定位  | Membrane,Single-pass type I membrane protein                                                  |
| 纯化    | Affinity purification                                                                         |

## 抗原信息


|      |                                                                                                                                      |
|------|--------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein corresponding to Human DCR2.                                                                              |
| 序列   | VLFRRRSCPSRVPGAEDNARNETLSNRYLQPTQVSEQEIQGQELAELTGVTVELPEEPQRLLEQAEAEGCQRRRLVPVND<br>ADSADISTLLDASATLEEGHAKETIQDQLVGSEKLFYEDEAGSATSCL |

## 靶点信息

|      |                                                                                                                                                                                                                                                                                                                                       |
|------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | The protein encoded by this gene is a member of the TNF-receptor superfamily. This receptor contains an extracellular TRAIL-binding domain, a transmembrane domain, and a truncated cytoplamic death domain . This receptor does not induce apoptosis, and has been shown to play an inhibitory role in TRAIL-induced cell apoptosis. |
| 基因ID | 8793                                                                                                                                                                                                                                                                                                                                  |

|       |                                              |
|-------|----------------------------------------------|
| 基因名   | TNFRSF10D                                    |
| Swiss | Q9UBN6                                       |
| 别名    | TNFRSF10D;CD264;DCR2;TRAIL-R4;TRAILR4;TRUNDD |

产品验证



Western blot analysis of DCR2 expressed in THP-1 using DCR2 Rabbit mAb at 1:1000. Secondary antibody : HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)