

Torsin A Rabbit mAb

货号: B29823

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB FC
推荐浓度	WB: 1:500 - 1:2000 FC: 1:20 - 1:50
理论分子量	22kDa/37kDa
实测分子量	35kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,BT-474,DU145
细胞定位	Cell projection,Cytoplasm,Cytoplasmic vesicle,Cytoplasmic vesicle membrane,Endoplasmic reticulum lumen,Nucleus membrane,Peripheral membrane protein,cytoskeleton,growth cone,secretory vesicle,synaptic vesicle
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Torsin A.
序列	KDLDNLFQGHLAKKIILNAVFGFINNPCKKPLTSLHGWTGTGKNFVSKIIAENIYEGGLNSDYVHLFVATLHFPHASNITLYKDQLQLWIRGNVSACARSIFIDEMDKMHAGLIDAIKPFLDYDLVDGVSYQKAMFIFLSNAGAERITDVALDFWRSGKQREDIKLKDIEHALSVSVFNKNKSGFWHSSLIDRNLIDYFVPFLPLEYKHLKMCIRVEMQSRGYEIDEDIVSRVAEEMTFFPK EERVFSKDGCKTVFTKLDYYYDD

靶点信息

研究背景	The protein encoded by this gene is a member of the AAA family of adenosine triphosphatases (ATPases), is related to the Clp protease/heat shock family and is expressed prominently in the substantia nigra pars compacta. Mutations in this gene result in the autosomal dominant disorder, torsion dystonia 1.
------	---

基因ID	1861
基因名	TOR1A
Swiss	O14656
别名	TOR1A;DQ2;DYT1;torsin-1A

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn