

Myb Rabbit mAb

货号: B29832

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	37-46kDa/63-85kDa
实测分子量	80kDa,72kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Myb.
序列	MARRPRHSIYSSDEDDDEFEMCDHDYDGLLPKSGKRHLGKTRWTRREEDEKLKKLVEQNGTDDWKVIANYLPNRTDVQC QHRWQKVLNPELIKGPWTKEED

靶点信息

研究背景	This gene encodes a protein with three HTH DNA-binding domains that functions as a transcription regulator. This protein plays an essential role in the regulation of hematopoiesis. This gene may be aberrantly expressed or rearranged or undergo translocation in leukemias and lymphomas, and is considered to be a n oncogene. Alternative splicing results in multiple transcript variants.
------	---

基因ID	4602
基因名	MYB
Swiss	P10242
别名	MYB;Cmyb;c-myb;c-myb_CDS;efg

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn