

Eg5 Rabbit mAb

货号: B29844

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IF/ICC IP
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	119kDa
实测分子量	119kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Jurkat,SW620
细胞定位	Cytoplasm,cytoskeleton,spindle pole
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Eg5.
序列	SVKHSDKLNNGNLEKISQETEQRCESLNTRTVYFSEQWVSSLNEREQELHNLLEVVSQCCEASSSDITEKSDGRKAAHEKQ HNIFLDQMTIDEDKLIAQNLELNETIKIGLTKLNCFLQDLKLDIPTGTTPQRKSYLYPSTLVRTEPREHLLDQLKRKQPELLM MLNCSENNKEETIPDVDVEEAVLGQYTEEPLSQEPSVDAGVDCSSIGGVPPFQHKKSHGKDKENRGINTLERSKVEETTEH LVTKSRLPLRAQINL

靶点信息

研究背景	This gene encodes a motor protein that belongs to the kinesin-like protein family. Members of this protein family are known to be involved in various kinds of spindle dynamics. The function of this gene product in cludes chromosome positioning, centrosome separation and establishing a bipolar spindle during cell mit osis.
------	---

基因ID	3832
基因名	KIF11
Swiss	P52732
别名	KIF11;EG5;HKSP;KNSL1;MCLMR;TRIP5

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn