

Caspase-2 Rabbit mAb

货号: B29855

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	12kDa/30kDa/48kDa
实测分子量	48kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,HeLa,293T
细胞定位	cytoplasm,cytosol,mitochondrion,nucleolus,nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Caspase-2.
序列	GKNHTQSPGCEESDAGKEELMKMRLPTRSDMICGYACLKGNAAMRNTKRGSWYIEALTQVFSEACDMHVADMLVKVN ALIKEREGYAPGTEFHRCKEMSEYCSLTCQLYLFPGYPPT

靶点信息

研究背景	This gene encodes a member of the cysteine-aspartic acid protease (caspase) family. Caspases mediate cellular apoptosis through the proteolytic cleavage of specific protein substrates. The encoded protein may function in stress-induced cell death pathways, cell cycle maintenance, and the suppression of tumorigenesis. Increased expression of this gene may play a role in neurodegenerative disorders including Alzheimer's disease, Huntington's disease and temporal lobe epilepsy. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jan 2011]
基因ID	835

基因名	CASP2
Swiss	P42575
别名	CASP-2; ICH1; NEDD-2; NEDD2; PPP1R57

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn