

ALKBH1 Rabbit mAb

货号: B29862

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	43kDa
实测分子量	43kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,293T,HepG2,MCF-7
细胞定位	Mitochondrion,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human ALKBH1.
序列	YQWHWVKQCLKLYSQKPNVCNLDKHMSKEETQDLWEQSKEFLRYKEATKRRPRSLLLEKLRWVTVGYHYNWDSKKYSA DHYTPFPSDLGFLSEQVAAACGFEDFRAEAGILNYRLDSTLGIHVDRSELDSKPLLSFSFGQSAIFLLGGLQRDEAPTA MFMHSGDIMMSGFSRLLNHAVPRVLPNPEGEGLPHCLEAPLPAVLPRDSMVEPCSMEDWQVCASYLKRTARVNMTVRQV LATDQNFPLEIEDEKRDISTEGFCHLDDQNSEVKRARINPDS

靶点信息

研究背景	This gene encodes a homolog to the E. coli alkB gene product. The E. coli alkB protein is part of the adaptive response mechanism of DNA alkylation damage repair. It is involved in damage reversal by oxidative demethylation of 1-methyladenine and 3-methylcytosine.
------	--

基因ID	8846
基因名	ALKBH1
Swiss	Q13686
别名	ALKBH1;ABH;ABH1;ALKBH;alkB;hABH

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn