

VAMP1 Rabbit mAb

货号: B29884

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	16kDa
实测分子量	16kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Rat brain
细胞定位	Cell junction,Cytoplasmic vesicle,Cytoplasmic vesicle membrane,Mitochondrion outer membrane,Single-pass type IV membrane protein,secretory vesicle,synapse,synaptic vesicle membrane,synaptosome,synaptosome
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human VAMP1.
序列	MSAPAQPPAEGTEGTAPGGGPPGPPNMTSNRRLQQTQAQVEEVVDIIRVNVDKVLERDQKLSELDDRADALQAGASQ FESSAAKLKRKYWWKNCKMMIMLGAICAIIVVVIVYFFT

靶点信息

研究背景	Synapotobrevins, syntaxins, and the synaptosomal-associated protein SNAP25 are the main components of a protein complex involved in the docking and/or fusion of synaptic vesicles with the presynaptic membrane. The protein encoded by this gene is a member of the vesicle-associated membrane protein (VAMP) /synaptobrevin family. Mutations in this gene are associated with autosomal dominant spastic ataxia 1. Multiple alternative splice variants have been described, but the full-length nature of some variants has not been defined. [provided by RefSeq, Jul 2014]
基因ID	6843
基因名	VAMP1
Swiss	P23763
别名	SPAX1; SYB1; VAMP-1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn