

MEF2A Rabbit mAb

货号: B29907

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB ICC
推荐浓度	WB: 1:500 - 1:2000 ICC: 1:50 - 1:200
理论分子量	45kDa/46kDa/53kDa/54kDa
实测分子量	55kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Raji,MCF7,Mouse spleen,C6
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human MEF2A.
序列	MGRKKIQITRIMDERNRQTLRKKGLNGCESPDADDYFEHSPLSEDRFSKLNEDSDFIKRGPPGLPPQNFSMSVTVPVTSP NALSYTNPSSLVSPSLAA

靶点信息

研究背景	The protein encoded by this gene is a DNA-binding transcription factor that activates many muscle-specific, growth factor-induced, and stress-induced genes. The encoded protein can act as a homodimer or as a heterodimer and is involved in several cellular processes, including muscle development, neuronal differentiation, cell growth control, and apoptosis. Defects in this gene could be a cause of autosomal dominant coronary artery disease 1 with myocardial infarction (ADCAD1). Several transcript variants encoding different isoforms have been found for this gene.
------	--

基因ID	4205
基因名	MEF2A
Swiss	Q02078
别名	MEF2A;ADCAD1;RSRFC4;RSRFC9;mef2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn