

HNF-4-alpha Rabbit mAb

货号: B29952

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	43-56kDa
实测分子量	53kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SW480
细胞定位	Nucleus
纯化	Affinity purification

抗原信息



抗原信息	Recombinant fusion protein corresponding to Human HNF-4-alpha.
序列	EWAKYIPAFCELPDDQVALLRAHAGEHLLLGATKRSMVFKDVLGNDYIVPRHCELAEMSRVSIRILDELVLPFQELQID DNEYAYLKAIIFDPAK

靶点信息

研究背景	The protein encoded by this gene is a nuclear transcription factor which binds DNA as a homodimer. The encoded protein controls the expression of several genes, including hepatocyte nuclear factor 1 alpha, a transcription factor which regulates the expression of several hepatic genes. This gene may play a role in development of the liver, kidney, and intestines. Mutations in this gene have been associated with monogenic autosomal dominant non-insulin-dependent diabetes mellitus type I. Alternative splicing of this gene results in multiple transcript variants encoding several different isoforms.
------	---

基因ID	3172
基因名	HNF4A
Swiss	P41235
别名	HNF4A;FRTS4;HNF4;HNF4a7;HNF4a8;HNF4a9;HNF4alpha;MODY;MODY1;NR2A1;NR2A21;TCF;TCF14

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn