

AGL Rabbit mAb

货号: B29961

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	172kDa/174kDa
实测分子量	175kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,SW620,Mouse heart,Mouse skeletal muscle
细胞定位	Cytoplasm
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human AGL.
序列	EGFNITAGVDEETGFVYGGNRFNCGTWMDKMGESDRARNRGIPATPRDGSAVEIVGLSKSAVRWLELSKKNIFPYHEVT VKRHGKAIKVSYDEWNRKIQDNFEKLFHVSEDPSDLNEKHPNLVHKRGIYKDSYGASSPWCDYQLRPNFTIAMVVAPELF TTEKAWKALEIAEKKLLGPLGPKTLDPDDMVYCGIYDNALDNDNYNLAKGFNYHQGPWLWPIGYFLRAKLYFSRLMGPE TTAKTIVLVKNVLSRHYVHLERSPWKGLPELTNENAQYCPFCSETQAWSIATILETLYDL

靶点信息

研究背景	This gene encodes the glycogen debrancher enzyme which is involved in glycogen degradation. This enzyme has two independent catalytic activities which occur at different sites on the protein: a 4-alpha-glucotransferase activity and a amylo-1,6-glucosidase activity. Mutations in this gene are associated with glycogen storage disease although a wide range of enzymatic and clinical variability occurs which may be due to tissue-specific alternative splicing. Alternatively spliced transcripts encoding different isoforms have been described.
基因ID	178
基因名	AGL
Swiss	P35573
别名	AGL;GDE

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn