

IRF3 Rabbit mAb

货号: B29969

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	12kDa/16kDa/33kDa/47kDa/49kDa
实测分子量	57kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,HeLa
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human IRF3.
序列	ASSLE

靶点信息

研究背景	This gene encodes a member of the interferon regulatory transcription factor (IRF) family. The encoded p rotein is found in an inactive cytoplasmic form that upon serine/threonine phosphorylation forms a compl ex with CREBBP. This complex translocates to the nucleus and activates the transcription of interferons al pha and beta, as well as other interferon-induced genes. Alternatively spliced transcript variants encodin g multiple isoforms have been observed for this gene.
------	---

基因ID	3661
基因名	IRF3
Swiss	Q14653
别名	IIAE7;IRF3

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn