

SRD5A2 Rabbit mAb

货号: B29976

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IP
推荐浓度	WB: 1:500 - 1:2000 IP: 1:20 - 1:50
理论分子量	28kDa
实测分子量	28kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	PC-3,Mouse liver,Mouse testis,Mouse kidney,Rat testis,Rat kidney
细胞定位	cell body fiber
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human SRD5A2.
序列	MQVQCQQSPVLG SATLV ALGALALYVAKPSGYGKHTE SLKPAATRLPARAAWFLQELPSFAVPAGILARQPLSLFGPPGT VLLGLFCVHYFHRTFVYSL

靶点信息

研究背景	This gene encodes a microsomal protein expressed at high levels in androgen-sensitive tissues such as the prostate. The encoded protein is active at acidic pH and is sensitive to the 4-azasteroid inhibitor finasteride. Deficiencies in this gene can result in male pseudohermaphroditism, specifically pseudovaginal perineoscrotal hypospadias (PPSH). [provided by RefSeq, Jul 2008]
基因ID	6716

基因名	SRD5A2
Swiss	P31213
别名	

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn