

KDM5A Rabbit mAb

货号: B29984

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IF/ICC IP FC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50 FC: 1:20 - 1:50
理论分子量	186kDa/192kDa
实测分子量	192kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,HeLa,MCF7
细胞定位	Nucleus,nucleolus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human KDM5A.
序列	

靶点信息

研究背景	This gene encodes a member of the Jumonji, AT-rich interactive domain 1 (JARID1) histone demethylase p protein family. The encoded protein plays a role in gene regulation through the histone code by specifically demethylating lysine 4 of histone H3. The encoded protein interacts with many other proteins, including r etinoblastoma protein, and is implicated in the transcriptional regulation of Hox genes and cytokines. This gene may play a role in tumor progression.
------	--

基因ID	5927
基因名	KDM5A
Swiss	P29375
别名	KDM5A;RBBP-2;RBBP2;RBP2

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn