

# ITK Rabbit mAb

货号: B29985

产品信息

|       |                                                                                               |
|-------|-----------------------------------------------------------------------------------------------|
| 反应    | Human                                                                                         |
| 宿主    | Rabbit                                                                                        |
| 克隆性   | Monoclonal                                                                                    |
| 预测反应  |                                                                                               |
| 应用    | WB IHC FC                                                                                     |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IHC:</b> 1:50 - 1:200<br><b>FC:</b> 1:20 - 1:50               |
| 理论分子量 | 71kDa                                                                                         |
| 实测分子量 | 72kDa                                                                                         |
| 形式    | Liquid                                                                                        |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                  |
| 阳性对照  | Jurkat,Molt-4,THP-1(Negative control)                                                         |
| 细胞定位  | Cytoplasm                                                                                     |
| 纯化    | Affinity purification                                                                         |

抗原信息


|      |                                                        |
|------|--------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein corresponding to Human ITK. |
| 序列   |                                                        |

靶点信息

|      |                                                                                                                                                                                                                                                  |
|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes an intracellular tyrosine kinase expressed in T-cells. The protein contains both SH2 and SH3 domains which are often found in intracellular kinases. It is thought to play a role in T-cell proliferation and differentiation. |
| 基因ID | 3702                                                                                                                                                                                                                                             |

|       |                          |
|-------|--------------------------|
| 基因名   | ITK                      |
| Swiss | Q08881                   |
| 别名    | ITK;EMT;LPFS1;LYK;PSCTK2 |

产品验证



Western blot analysis of ITK expressed in Jurkat,Molt-4,THP-1(Negative control) using ITK Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)