

Phospho-Paxillin (Y118) Rabbit mAb

货号: B30040

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	68kDa
实测分子量	68kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-431,NIH/3T3,C6
细胞定位	Cell junction,Cytoplasm,cell cortex,cytoskeleton,focal adhesion
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Phospho-Paxillin (Y118).
序列	HVYSF

靶点信息

研究背景	This gene encodes a cytoskeletal protein involved in actin-membrane attachment at sites of cell adhesion to the extracellular matrix (focal adhesion). Alternatively spliced transcript variants encoding different iso forms have been described for this gene. These isoforms exhibit different expression pattern, and have di fferent biochemical, as well as physiological properties (PMID:9054445). [provided by RefSeq, Aug 2011]
基因ID	5829
基因名	PXN

Swiss	P49023
别名	PXN; paxillin

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn