

ELK1 Rabbit mAb

货号: B30056

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC ICC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 ICC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	11kDa/44kDa
实测分子量	62kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Mouse placenta,Rat heart
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human ELK1.
序列	PSPLEACLEAEEAGLPLQVILTPPEAPNLKSEELNVEPGLGRALPPEVKVEGPKEELEVAGERGFVPETTKAEPEVPPQEGVPA RLPAVVMDTAGQA

靶点信息

研究背景	This gene is a member of the Ets family of transcription factors and of the ternary complex factor (TCF) s ubfamily. Proteins of the TCF subfamily form a ternary complex by binding to the the serum response fact or and the serum response element in the promoter of the c-fos proto-oncogene. The protein encoded by this gene is a nuclear target for the ras-raf-MAPK signaling cascade. This gene produces multiple isoforms by using alternative translational start codons and by alternative splicing. Related pseudogenes have bee n identified on chromosomes 7 and 14.
------	--

基因ID	2002
基因名	ELK1
Swiss	P19419
别名	ELK1

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn